



Scan the QR code for more information

Multi-Drug Screen Urine Test for Prescription Use Only



● Suitable For

- ✓ Medication monitoring
- ✓ Pain management applications etc
- ✓ Criminal justice
- ✓ Point of care sites

Drug	Cut Off (ng/mL)
MET	500/1000
MOP/OPI	300/2000
MTD	300
OXY	100
PCP	25
PPX	300
TCA	1000
THC	20/50
6-MAM	10

Drug	Cut Off (ng/mL)
AMP	500/1000
BAR	300
BUP	10
BZO	300
COC	150/300
EDDP	300
MDMA	500
FEN1 *	1
FEN5 *	5

* Note: Fentanyl 1 and Fentanyl 5 are currently in the registration process and will be available soon.

● Product Format: Cup / Dip Card

10140 Mesa Rim Rd, Unit 1, San Diego, CA 92121

www.advinbio.com · E-MAIL:support@advinbio.com · TELEPHONE:1-858-866-8382