




Scan the QR code for more information

# Vinscreen Oral Tests



## Saliva Tracer

- ★ Minimal Saliva Needed
- ★ Rapid Sample Collection
- ★ Detects Parent THC
- ★ Simple and Intuitive Design



## O-Flow

- ★ Detect Multiple Drugs
- ★ Simple One-Step Process
- ★ Lightweight and Portable
- ★ Alcohol Testing Available

Drug	Cut Off (ng/mL)	Drug	Cut Off (ng/mL)
AMP	25/50	MQL	100
BAR	50	OPI/MOP	10/30/40/50
BUP	5	OXY	20/40
BZO	10/20	PCP	10
COC	10/15/20/50	PPX	30/50
COT	30/50	SMA/K2	25
FYL	10	SMP	10
KET	30/50	THC	15/40/50
6-MAM	3/5/10	TML	30/50
MDMA	50	ZOP	20
MET	25/50	ALC	0.02%(20mg/dL)
MTD	30		

10140 Mesa Rim Rd, Unit 1, San Diego, CA 92121

www.advinbio.com · E-MAIL:support@advinbio.com · TELEPHONE:1-858-866-8382